SAFEGUARDMEDICAL™

Prometheus Traction Splint

Contents: 1Pcs

Safeguard Medical Technologies Ltd.
The Old Rectory
Hope Under Dinmore
Herefordshire, HR6 0PW
United Kingdom

DTF050_0003_Rev7_05/21

Authorized Representative
in the EC

Safeguard Technologies Ltd
Willow Grove, Delgany, A63 XY89
Co Wicklow, Ireland

Prometheus Traction Splint

Intended Use

The Prometheus Traction Splint is indicated for use on suspected limb fractures.

User Profile

The Prometheus Traction Splint is designed to be used by medical professionals in the pre-hospital setting for immobilisation of suspected limb fractures.

Patient Profile


The Prometheus Traction Splint is designed to be used on those patients with a suspected long bone fracture. Proper use can reduce haemorrhage, decrease pain and associated fracture complications. The femur is the longest and strongest bone in the human body; the surrounding muscles are very strong. When the femur fractures, the surrounding muscles can spasm, resulting in fracture displacement, which causes severe pain, haemorrhage (internal or external), muscle and neurovascular injury. The intended therapeutic indication of a traction splint is to reduce haemorrhage, relieve pain and to stabilize the fracture site, preventing further injury. The traction splint can be in place for up to 24 hours.

Performance Characteristics

- The pole can be adjusted for length
- Robust and durable
- Straps can be extended to cater for an ever increasingly sized population
- Safe
- Capable of required traction
- Does not impede access to groin or pelvis
- Requires only single operator
- Not ferromagnetic
- Lightweight

Contents

- AH - Ankle Hitch
- TL - Thigh Loop
- P - Carbon Fibre Pole (Lower Thigh Strap pre-fitted to pole)
- LLS - Lower Leg Strap
- UTS - Upper Thigh Strap

Clinical Benefits

Reduction of a displaced fracture reduces haemorrhage into the fracture site. Patients usually experience significant relief from pain after the limb is restored to normal alignment.

Contraindications

Do not use on neck of femur (i.e. hip), ankle or knee injuries.

Residual Risks/Warnings

- Hold the poles away from the joints being careful not to catch fingers as it can cause injury.
- Do not overtighten any strap on the device in order not to restrict blood flow or cause skin breakdown/nerve damage. Monitor the distal pulse regularly to detect any restricted blood flow and loosen straps if required.
- If there is any sign of patient discomfort from the Thigh Loop, apply loose padding under the strap against the inner thigh.
- The product is intended for single use only. Re-use will risk cross infection and loss of product functionality.
- After use, the device should be disposed of as clinical waste as per local protocols, failure to do so will result in cross contamination.
- Do not use after expiry date. Use after expiry date will risk loss of product functionality.
- Any serious incident that has occurred in relation to this device should be reported to the manufacturer (vigilance@safeguardmedical.com) and applicable national competent authority.


Instructions for Use

1. Apply the ankle hitch around the ankle with the stirrup strap beneath the foot

2. Secure the ankle hitch with the 'hook and loop' strap.


3. Tighten the stirrup by pulling on the green tab.

4. Apply the Thigh Loop by unclipping the buckle, feeding the male part of the buckle behind the knee and slide up the leg.

5. Clip the buckle together and tighten the strap so that the pole receptacle is positioned approximately at the patient's belt line. The strap should not be so tight as to restrict blood flow. Thread the strap so the buckle is on the anterior thigh.

6. Snap out the pole carefully to ensure that each joint is securely in place. **NOTE:** Hold the poles away from the joints being careful not to catch fingers.

7. Lay the pole alongside the patient's leg with the dart at the foot end. Approximately one section of this pole should extend below the foot.

8. The top end of the pole should now be sized against the receptacle. Unsnap sections of the pole as required and fold over.


9. Firmly insert both sections into the receptacle.

10. Loosely apply the Lower Thigh Strap by feeding the strap under the knee and fastening with the 'hook and loop' material.

11. Fit the yellow tab over the dart at the foot end of the pole.

12. Stabilise end of pole by resting against knee. Apply traction to the patient's lower leg. Once the appropriate amount of traction has been applied, tension the traction strap by pulling on the red tab. The degree of tension should be in accordance with local protocols.

13. Apply the Upper Thigh Strap by placing it underneath the knee and sliding up towards the top of the thigh. Secure the strap using the 'hook and loop' material.

14. Adjust the Lower Thigh Strap as necessary. (Ideally this should rest just above the patella).

15. Place the Lower Leg Strap between the ankle hitch and mid-calf and fasten using the 'hook and loop' material.

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