


Team-8® Stretcher



Overview

The Team-8® Stretcher is a lightweight, compact, rapidly deployable stretcher with colour coded webbing straps and quick release buckles to get casualties out of a hazardous or hostile environment as quickly and safely as possible. Made of lightweight, durable polyethylene, it is designed to be used in urban and non-urban environments. The Team-8® Stretcher is able to be deployed rapidly from its carrying bag, and can be used as a lift and a drag stretcher.

Key Benefits

-Lightweight & compact, for quick casualty evacuation.

- Safe.
- Lifiable.
- Robust.
- Product code: PRO188NC

Specifications



Dimensions: 210.00cm x 46.50cm x 0.10cm (Pack dimensions: 56.00cm x 16.00cm)



Weight: 2.40kg



Product shelf life: 5 years