

# 68 MEDICAL SOLUTIONS

2026 PRODUCT GUIDE  
[WWW.68MEDICAL.COM](http://WWW.68MEDICAL.COM)

# VETERAN OWNED VETERAN OPERATED



## URIAH POPP

FOUNDER, CEO

US Army Special Operation and  
Special Mission Unit Medic  
Inventor/ Product Developer



## RYLAN MILLER

PARTNER, CFO

US Navy Special Warfare  
Tactical Surveillance Element Operator  
CEO Zenith investment group



## SHINTELL POPP

PARTNER, VP QARA


US Army Transportation Specialist  
Cadaveric Lab Manager



## PAUL STULTZ


PARTNER, COO

US Army Cannon Crewmember  
Field Artillery Targeting Technician  
GOV Acquisition Level 3



6:8 MEDICAL SOLUTIONS

# PROUD DISTRIBUTOR FOR



## Mission:

To revolutionize prehospital emergency care by delivering innovative, life-saving medical devices that address critical gaps in care, enhance operational efficiency, and improve patient outcomes in both civilian and military settings.

## Vision:

To be the global leader in prehospital medical technology, pioneering transformative solutions that empower healthcare providers, save lives, and set new standards in trauma care.

## Purpose:

To advance prehospital medicine through the development of cutting-edge devices and targeted solutions for managing life-threatening conditions, ensuring that no patient suffers from preventable injury or death.

# OUR GLOBAL



**PROMOTED**

artzena



CSD






BLACKSMITHS

DISTRIBUTION PARTNERS

6:8 MEDICAL SOLUTIONS 8:9


# POPP POINT PRESSURE DEVICE P3D TECHNOLOGY



## Point Pressure versus Constriction

- Provide perfusion from non-targeted vessels
- Mitigate painful venous engorgement
- Deliberate devices for tourniquet conversion
- Limit distal ischemia
- Avoid compartment syndrome

LIMITS IMPROVISING,  
REBLEED,  
AND  
UNNECESSARY BLOOD LOSS



6:8 MEDICAL SOLUTIONS



# BLEEDING CONTROL

6:8 MEDICAL SOLUTIONS

## XMETIX TAK 710 Tourniquet

PN: TAK 710

FDA/ CE

Single Patient Use

Latex-Free

5-year shelf life

Dimensions:

37" L X 1.5" W



## CAT Tourniquet

PN: 30-0001

FDA/ CE

CoTCCC Recommended

Single Patient Use

Latex-Free

5-year shelf life

Dimensions:

37" L X 1.5" W



TMT

PN: 60064

CoTCCC Recommended

FDA

Single Patient Use

Latex-Free

10-year shelf life

Dimensions:

38.7" L X 2.25" W



## READY WRAP

PN: RDYWRP

FDA /CE

Single Patient Use

High Stretch **Grip8** Rubber Wrap

Latex-Free and Powder Free

10-year shelf life

Dimensions:

4" X 50"



**PERFECT FOR K9 USE**

6:8 MEDICAL SOLUTIONS

# HEMGUARD

WORLD'S **FIRST** TARGETED PRESSURE TRAUMA BANDAGE




Your Life May Depend,  
On What You Packed



6:8 MEDICAL SOLUTIONS

# HEMGUARD

## TARGETED PRESSURE



**PN: HGTPD-04T**

**NSN: 6515-01-732-6982**

**HGTPD-06T**

**NSN: 6515-01-732-9490**

**FDA/ Trade Compliant**

**PN: HGTPD-04 / HGTPD-06**

**FDA/ CE**

**Single Patient Use**

**Latex-Free**

**5-year shelf life**


**Packaged Dimensions:**

HGTPD-04 4.5" x 9.75" x 1.5"


HGTPD-06 4.5" x 10.75" x 1.5"

**Weight: 4 oz**

## EMERGENCY CONSTRICTION



**REDIRECTION CHANNELS SIMPLIFY BANDAGING  
WHILE RISER TARGETS PRESSURE AND PREVENTS REBLEED**



# HEMGUARD<sup>®</sup> FLAT

PN: HGFTPD-4 / HGFTPD-6

FDA/ CE

Single Patient Use

Latex-Free

5-year shelf life

Packaged Dimensions:

HGTPD-F4 4.5" x 8.75" x 1"

HGTPD-F6 6.5" x 9.75" x 1"

Weight: 4 oz



**4 Way Redirection;** Anchor Bandage to Maximize Pressure

**Mechanical Advantage;** In 2 Directions if Desired



NOW IN A **LOW-PROFILE** DESIGN FOR COMPACT STORAGE



**6:8 MEDICAL SOLUTIONS**

# THE ADVANTAGE

PRESSURE BAR AND PRESSURE BAR SET

SUPERFICIAL



MID LAYER



DEEP



TARGETED CONTROL



WORKS WITH  
WHAT YOU ARE  
ALREADY  
CARRYING



# Pressure Bar Set

PN: TADV-01

NSN: 6515-01-733-0896

# Pressure Bar

PN: TADV-02

NSN: 6515-01-733-6096



## WHY TARGETED PRESSURE?

Maintains Constant Pressure Prevents Rebleed and Supports Clot Formation

Allows For Distal Perfusion and Venous Return Preventing Constriction Injury and

**HOLDING PRESSURE**  
**SO, YOU DON'T HAVE TO**

Single Patient Use  
8 YR Shelf life  
Latex-Free  
CE/ FDA  
Packaged Dimensions  
5.0" x 3.75" x 1.5"  
4" x 2" x 1.5"  
Weight  
3 oz / 5 oz  
Made in the USA



EVER TRIED TO OPEN A PACKAGE  
WITH THESE HANDS?

**PN: CZFG-01**

**FDA/ CE**

**Single Patient Use**

**Latex-Free**

**5-year shelf life**

**Dimensions:**

4.5" X 4.1 yards

11.5 cm X 3.75 meters

**Weight: 3 oz**

# COMPACT Z-FOLD GAUZE



**Packaged with the  
Professional in Mind**

6:8 MEDICAL SOLUTIONS

# SEALSKIN

WOUND WRAP

CLEAR


EASY TEAR

NON-ADHERENT

NON-ADHESIVE

WATERPROOFING

SEALING WRAP



**PN: SEALSKN**

FDA

Single Patient Use

Latex-Free

5-year shelf life

Dimensions: 3" X 50'

Weight: 3 oz

Made in the USA



# CHITOGAUZE

**XR PRO**

COMPRESSED Z-FOLDED CHITOSAN BASED HEMOSTATIC GAUZE



**PN: 1090**

CoTCCC Recommended

FDA/ CE

Single Patient Use

Latex-Free

5-year shelf life

3" X 4 YRDS

Professionals who work in these environments,  
Demand products that will work there.


PLANT BASED SAFE AND EFFECTIVE

ankaferd  
**BLOODSTOPPER**<sup>®</sup>



Creates a specific erythrocyte protein network that physiologically directs hemostasis in plasma and serum in a short time.


CE **External** and **Internal** Use



**250** PUBLISHED ARTICLES AND STUDIES IN INTERNATIONALLY ACCEPTED PUBLICATIONS.

OVER **14** YEARS OF USE

**FINALLY.**  
A Solution for Both Compressible and Non-Compressible Hemorrhage



**6:8 MEDICAL SOLUTIONS**

NOT AVAILABLE IN EVERY COUNTRY

[1] Alpay A, Evren C, Bektaş S, Uğurbas Canturk S, Uğurbay SH, Cinar F. Effects of the folk medicinal plant extract Ankaferd Blood Stopper® on the ocular surface. Cutaneous and Ocular Toxicology, 2011; 30(4): 280-285. Cutaneous and Ocular Toxicology, 2011; 30(4): 280-285




# AIRWAY AND RESPIRATORY MANAGEMENT

6:8 MEDICAL SOLUTIONS

# PROCRIC

CRICOTHYROIDOTOMY KIT

6.0



Patented Rota-Trach  
Rotational Flange;  
Prevents dislodgement  
while ventilating the  
patient



6.0mm Tube  
Obturator  
Neck Strap  
Syringe  
Scalpel  
Wedge

PN: PROCRIC

NSN 6515-01-737-1154

Single Patient Use

Latex-Free

3-year shelf life

Packaged Dimensions:

9.0" x 6.5" x 1.5"

Weight: 3 oz

CE/ FDA

TAA Compliant

Patented **Popp Obturator**; Integrates a rounded flat tipped introducer to **Ease Insertion** through the incision and a syringe holder to **Streamline Procedure**



6:8 MEDICAL SOLUTIONS

# SURSECUR

™ 28 FR NPA

NASOPHARYNGEAL AIRWAY



NO TAPE NEEDED

PATENT PENDING

PN: SURNPA / SURNPA-T

NSN 6515-01-737-1154

CE/ FDA

Single Patient Use

Latex-Free

5-year shelf life

Packaged Dimensions:

8.5" x 2.75" x .5"

Weight: 3 oz

TAA Compliant

## Anti-Rotation SurSecur Clip

Prevents rotating, backing out, and indicates bevel orientation

## Posterior Fenestrations

Allows for passing air in the event the lumen becomes blocked

Simplifying  
Airway

Management



**PN: PDF0680**

**Single Patient Use**

**Latex-Free**

**5-year shelf life**

**Packaged Dimensions:**

**5.0" x 6.0"**

**FDA**

**TAA Compliant**

**PN: DCN-01**

**Single Patient Use**

**Latex-Free**

**5-year shelf life**

**Packaged Dimensions**

**:**

**6.0" x .5"**

**FDA/CE**



# HYPOTHERMIA MANAGEMENT

A person is lying on their back on a snowy, icy beach. They are wearing a bright orange survival suit with black straps. The person's head is covered with a black hat. The beach is covered in snow and ice, with the ocean in the background. The sky is a mix of blue and orange, suggesting a sunset or sunrise. The overall scene is cold and desolate.

6:8 MEDICAL SOLUTIONS

# LITE FIGHTER

POLELESS LITTER

75G HEAT REFLECTIVE INSULATED CENTER PAD

RIPSTOP NYLON

3 G BUCKLE SECURING STRAPS

8 CARRY HANDLES

HOOK AND LOOP CENTER OPEN

18" W x 72" L



12" x 5"

Plate Carrier Hanger or  
Molle Attachment  
Storage Roll

PN: LFPTL

6:8 MEDICAL SOLUTIONS

# HEAT

## HYPOTHERMIA EMERGENCY ASSESSMENT AND TRANSPORT



**Flat Fold Litter**  
19" Wide 77" Long  
4 Securing Straps  
8 Carrying Handles

Head elevating Drag Hole

HEAT COMPLETE 3-PART SYSTEM  
FOLDABLE DRAG LITTER  
CLEARVIEW PATIENT CARRIER  
STORAGE BACKPACK

PN: HEAT  
PN: FLDS  
PN: CVPT  
PN: STRPK

250G ARCTIC  
150G TEMPERATE



### Storage Backpack

Flat Design  
12" Width  
19.5" Length  
3" Wide  
Laser Cut  
Molle Front

AVAILABLE IN:  
MULTICAM GREEN  
OLIVE DRAB GREEN  
RESCUE ORANGE



6:8 MEDICAL SOLUTIONS

# A SYSTEM FOR SIMULTANEOUS TREATMENT, TRANSPORT, AND HEAT LOSS PREVENTION



Upper Extremity Access



ClearView Chest Panel



Inguinal Access



Quilted Heat Reflective Lining




Lower Extremity Access



Dual Full Length Zipper Closure





MADE IN ITALY



# STILA TWIN

The rapid disconnect mechanism and premium carbon fiber composition of the Stila Twin allow the stretcher to drastically decrease its size and weight while maintaining superior strength and required capabilities.

The conformation and folding method of Stila Twin allow all accessories to be contained within it, making it “ready to use”.

Made in Italy



## Specifications:

Length: 186 cm

Width: 44 cm

Thickness: 0,4 cm

Weight: 4,6 kg

Folded length: 100 cm

Folded thickness: 9 cm

Weight: table 6.5 kg ,

with accessories 7.4 kg

Load capacity (EN 1865-

1:2015): 180 Kg



# FLUITANS

Fluitans allows the use of N-E-XT board as a rescue extractor in water environments, allowing the rescue and recovery of endangered or wounded people with rope or winch maneuvers (together with Tyrah rescue triangle).



Specifications:

Weight : 97 g without cylinder

33 g CO2 cylinder weight = 136 g

Length : 70 cm

Width : 37 cm

Maximum thickness (folded) = 3,6 cm

Buoyancy : 150N

# SAR-TFA


SAR-TFA is the floating aid kit to be used in combination with the TYRAH AR and the N-E-XT for the recovery of an injured patient in the aquatic environment also by lifting (Tyrah AR). The SAR-TFA was created to help the floating and recovery of patients from aquatic environments, guaranteeing all the handling precautions typically applied to traumatized patients.

Specifications:

Buoyancy aid : 50N

Weight : 700g

Width : 120cm



# N-E-XT

Emergency spine immobilization device that can be used in extrication or in combination with human external cargo harnesses for hoist or rope operations.

Specifications:

Length: 35,4 in (90cm)

Width: 11,8 in (30 cm)

Thickness: 0,16 in (0,4 cm)

Weight: 3,3 lb (1,5 kg)

# TYRAH

The Tyrah guarantees the immobilization and the support of the spine during rope maneuvers using a hoist.

The dedicated recovery system ensures the best position of the patient during handling and avoids problems related to compartment syndromes.



Specifications:  
Length : 22 cm  
Width : 22 cm  
Weight : 0,25 Kg

## N-E-X-T BAG

The NEXT system with its accessories can be transported with the dedicated hoistable bag or in standard soft bag.

**NORTHWALL®**  
INNOVATION




Specifications:  
Weight : 760 g  
Width : 114 cm  
Height : 116 cm  
Working load : 20 Kg to 150 Kg

## THE QUICK-HEAD

The Quick-Head immobilizer system is able to limit adjustments by adapting to a wider range of patients.

The Quick-Head immobilizer system with pre-calibrated elastics allows for safe operation with the assurance of not generating harmful compression and movement.



VI VIII  
**6:8 MEDICAL SOLUTIONS**

# ARCTIC WOLF

HAND WARMER



ClearView Window  
Protective Cover  
Soft Shell Exterior  
Zippered Expander



Heat Reflective  
Quilted Lining  
Rear Zipper Pocket  
Plate Carrier Hangers  
and Neck Strap

**PN: AWHW**

The ClearView Window allows you to use fine motor skills with your hands and provides you a field of view free from the cold.



**6:8 MEDICAL SOLUTIONS**



# ARCTIC WOLF XL

PN: AWHWX

Protective Removable Cover  
ClearView Window  
Cordura and Soft-Shell Exterior

Integrated Neck Strap  
Front Zipper Opening

Mesh Inner Pocket  
for Warming Pad  
Heat Reflective  
Quilted Lining

\*Perfect for Tablets and Larger Controllers

# 6:8 MEDICAL SOLUTIONS

# HEMOVAULT

BLOOD TRANSPORT SYSTEM



TEMPAID  
Cold Chain Packaging

ISTA 7D tested against Summer & Winter profiles.  
72 HRS+ of 2C- 8C Temperature Stability


MADE IN CANADA




2°C - 8°C  
15°C - 25°C  
-15°C to -25°C

 3L & 5L

 Long Duration  
Up to 72+ Hours

 Portable &  
Easy to Carry



PCM Plates fit both the 2-unit and 4-unit systems



- VIP LID NESTED INTO TRANSPORT BAG
- 6 PHASE CHANGE MANAGEMENT (PCM) RIGID and SEALED PLATES
- PATENTED ROTECTIVE CLEANABLE REINFORCED INNER TUB
- VIP PANELS NOW WITH REINFORCED CORNERS
- HEAVY DUTY 1000D EASY TRANSPORT BAG WITH CARRYING STRAP AND DUAL BUCKLE AND ZIPPER CLOSURE

# SUBCOMPONENT REPLACEABLE CONSTRUCTION

PN: HEMVLT-EMS4  
PN: HEMVLT-MIL4

PN: HEMVLT-EMS2  
PN: HEMVLT-MIL2

PN: QIF03-BUA1000



HemoTemp® II

1-3 4-6 7-9

CELSIUS

Smart Temperature Indicator Strips for Placement on Blood Bags

10°C

spotSee

# HYPOTHERMIA PREVENTION KIT



Includes:  
8 Cell HAWK Warming Grid  
3 Layer Blizzard Blanket

**PN: HYPO-BK3F-O**

[WWW.68MEDICAL.COM](http://WWW.68MEDICAL.COM)

# MEDICAL RESUPPLY KITS

## Your Choice of Tourniquet

CAT/TMT/TAK710 Tourniquet

CoTCCC Recommended Items  
Covering the **PREVENTABLE** causes  
of death from trauma.

## Key Products Included:

### Basic Kits Include:

Hemguard 4" Bandage (SurGrip Packaged)

(2) Compressed Z-Fold (SurGrip Packaged)

The Advantage Pressure Bar (SurGrip Packaged)

Sursecur NPA 28F (SurGrip Packaged)

Russell Chest Seal

Pressure Wrap

Mylar Blanket

Marker

Shears

### Advanced Kits Include:

Chitogauze XR PRO Hemostatic gauze

14G 3.25" Decompression Needle

TCCC Card

### Cold Weather Kits Include:

Blizzard 2-layer Blanket PBS2-370

Each product was selected after years of clinical and training evaluation to ensure efficacy when you need it most.

\*KITS COMES VACUUM SEALED



# MEDICAL RESUPPLY KITS

## IFAK Resupply Advanced CAT PN: IFAKA-CAT



### Kit Include:

- CAT Tourniquet
- Hemguard 4" Bandage
- Chito gauze XR Pro
- 14G Decompression Needle
- Compressed Z-Fold
- The Advantage Pressure Bar
- Russell Chest Seal
- SurSecur NPA 28F Marker
- Pressure Wrap
- Trauma Shears
- Mylar Blanket
- TCCC Card

## IFAK Resupply Basic CAT PN: IFAKB-CAT

### Kit Include:

- CAT Tourniquet
- Hemguard 4" Bandage
- (2) Compressed Z-Fold
- The Advantage Pressure Bar
- Russell Chest Seal
- SurSecur NPA 28F
- Pressure Wrap
- Mylar Blanket
- Marker
- Trauma Shears



\*KITS COMES VACUUM SEALED

# MEDICAL RESUPPLY KITS

## IFAK Resupply Advanced TAK PN: IFAKA-TAK

### Kit Includes:

- TAK 710 Tourniquet
- Hemguard 4" Bandage
- Chito gauze XR Pro
- 14G Decompression Needle
- Compressed Z-Fold
- The Advantage Pressure Bar
- Russell Chest Seal
- SurSecur NPA 28F Marker
- Pressure Wrap
- Trauma Shears
- Mylar Blanket
- TCCC Card



## IFAK Resupply Basic TAK PN: IFAKB-TAK

### Kit Include:

- TAK 710 Tourniquet
- Hemguard 4" Bandage
- (2) Compressed Z-Fold
- The Advantage Pressure Bar
- Russell Chest Seal
- SurSecur NPA 28F
- Pressure Wrap
- Mylar Blanket
- Marker
- Trauma Shears



\*KITS COMES VACUUM SEALED

# MEDICAL RESUPPLY KITS

## IFAK Resupply Advanced TMT PN: IFAKA-TMT



### Kit Include:

- TMT Tourniquet
- Hemguard 4" Bandage
- Chito gauze XR Pro
- 14G Decompression Needle
- Compressed Z-Fold
- The Advantage Pressure Bar
- Russell Chest Seal
- SurSecur NPA 28F
- Mylar Blanket
- Pressure Wrap
- Marker
- TCCC Card
- Trauma Shears

## IFAK Resupply Basic TMT PN: IFAKB-TMT

### Kit Include:

- TMT Tourniquet
- Hemguard 4" Bandage
- (2) Compressed Z-Fold
- The Advantage Pressure Bar
- Russell Chest Seal
- SurSecur NPA 28F
- Pressure Wrap
- Mylar Blanket
- Marker
- Trauma Shears



\*KITS COMES VACUUM SEALED

# MEDICAL RESUPPLY KITS

## Cold Weather IFAK Resupply Advanced CAT PN: IFAK-CWRS-CAT



### Kit Include:

- TMT Tourniquet
- Hemguard 4" Bandage
- Chito gauze XR Pro
- 14G Decompression Needle
- Compressed Z-Fold
- The Advantage Pressure Bar
- Russell Chest Seal
- Blizzard 2 Layer
- SurSecur NPA 28F
- Trauma Shears
- Pressure Wrap
- TCCC Card

## Cold Weather IFAK Resupply Advanced TMT PN: IFAK-CWRS-TMT

### Kit Include:

- TMT Tourniquet
- Hemguard 4" Bandage
- (2) Compressed Z-Fold
- The Advantage Pressure Bar
- Russell Chest Seal
- SurSecur NPA 28F
- Pressure Wrap
- Blizzard 2 Layer
- Marker
- Trauma Shears



\*KITS COMES VACUUM SEALED

# MEDICAL RESUPPLY KITS

## MICRO BLEEDER KIT PN: MCBLDR-BLK

Includes:

- TMT
- Advantage Pressure Bar
- Compressed Gauze
- RAP TQ Pouch



## DBAG IFAK PN: IFAK-ADV-TMT

Color Choice:  
 Multicam Green  
 Multicam Black  
 Multicam White  
 Black  
 Green  
 Tan  
 Red



### Kit Include:

- TMT Tourniquet
- Hemguard 4" Bandage
- Chito gauze XR Pro
- 14G Decompression Needle
- Compressed Z-Fold
- Russell Chest Seal
- SurSecur NPA 28F
- Pressure Wrap
- Mylar Blanket
- Trauma Shears
- Marker
- Tape
- TacLight



# M11

## MEDICAL BAG

### FEATURES

PN: M11MB-MC/ RO/ ODG

Laser Cut Molle

4 Exterior Tourniquet Holders

3 Removeable Inner Clear Pouches

Outer Zipper Pocket

Cordura Material

Rear Foldable Drag Litter Pocket



# 6:8 MEDICAL SOLUTIONS





Rear Flat Fold Litter Pocket

PATENTED QUAD KEY LOCK SYSTEM

COLOR OPTIONS AVAILABLE

## DIMENSIONS

Exterior

20" L X 12" W X 4.5" H

3 Inner Removable Pouches

11" L X 2.5" W X 5" H

3 Inner Mesh Pouches

11" L X .5" W X 6" H

Lower Hanger Pouch

11" L X 3" W X 4" H

# THE ORGANIZER

RANGE OR GO BAG,

AUTO OR ATV,

DESK OR DRAWER

## PN: ORGZR FEATURES

- 8 Elastic Pouches 4 Per Side
- Large Carrying Handle
- 2 Hook and Loop Patch Panels

6:8 MEDICAL SOLUTIONS 8:9

# DBAG IFAK

## FEATURES

### PN: DBAG-Color Code

- Molle Quick Detach
- Laser Cut Inner
- Dual TQ Side Elastic
- Two Pockets Under Laser Inner
- Diagonal Outer Pocket
- Cordura Material
- Rubber handle Pulls



### AVAILABLE:

Colors: Multicam / Multicam BLK/  
Multicam WHT/ Red/ OD GRN/  
Coyote Tan/ Black

6:8 MEDICAL SOLUTIONS



# 3MAG

VAMPIRE POUCH



BLOOD KIT STORAGE SOLUTION

MULTICAM / COYOTE TAN/ OLIVE DRAB / BLACK

6:8 MEDICAL SOLUTIONS

# RAPID ACCESS POUCH SIDE PULL

## FEATURES

PN: RAP-SP-MC

- Molle and Belt Attachment Rear
- Laser Cut Molle Front
- Large Side Pull Handles
- Hook and Loop Closure
- Inner Pouch Organizer

L X W X H



AVAILABLE:

Colors:  
Multicam  
Black  
Green

# 6:8 MEDICAL SOLUTIONS

# RAPID ACCESS POUCH STRAP HANGER



## FEATURES

PN: RAP-SH-MC

- 5.0" x 8" x 3" Outer Dimensions
- Hook and Loop 3" x 5" hanger panel
- Laser cut molle front with zipper pouch
- Lower tourniquet retention loop
- 90 degree open
- 2 staggered upper pouches with scissor holder
- 1 lower under pouch with elastic retention loops above



# 6:8 MEDICAL SOLUTIONS



BACK



FRONT



# SURGRIP PACKAGING



SURGRIP INDEXING HOLES

LOW LIGHT INDICATORS RINGS

CENTER TEAR ASSIST

PRODUCT VIEW WINDOW



**NVG COMPATIBLE**

**INNOVATION**  
WITH ATTENTION TO DETAIL

SWAY

## GIVE US A FOLLOW

six8\_medical\_solutio...    

Note...



**230** posts   **8,206** followers   **1,634** following

**6:8 Medical Solutions | FDA Medical Device MFG**  
Product/service  
Inventor of P3D Targeted Pressure Devices| ProCric|  
SurSecur NPA| Henguard Bandages| HEAT  
Hypothermia . Former SOF/SMU Medic... more  
[68medical.com](http://68medical.com)

**Professional dashboard**  
721.0K views in the last 30 days.

[Edit profile](#)   [Share profile](#)   [Email](#)

